

File D:\DATA

Injection Date : 28 Feb 15 2:45

Sample Name : 96206-92-7

Acq. Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Analysis Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Sample Info : *** No target masses specified ***

Method Info : Mobile Phase: A: water(10mM Ammonium hydrogen carbonate) B: ACN





Gradient: 5%-95% B in 1.5min

Flow Rate: 2.0ml/min

Column: XBridge C18, 4.6*50mm, 3.5um

Oven Temperature: 40 C

MASS Range:110-800



Integration Results for ELS1 A, Voltage

RetTim	Width	Area	Height	Area%
--------	-------	------	--------	-------

Integration Results for DAD1 A, Sig=214,4 Ref=off

RetTim	Width	Area	Height	Area%
2.11	0.02	3409.04	2674.93	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
2.11	0.02	3456.92	3161.06	100.00

Ret. Time: 2.10


<<<< POSITIVE SPECTRA >>>>




Ret. Time: 2.25




Ret. Time: 2.34




Ret. Time: 2.47




Ret. Time: 2.51



Ret. Time: 2.60



Ret. Time: 2.64



Ret. Time: 2.70

