

File D:\DATA\

Acq. Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Analysis Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Sample Info : *** No target masses specified ***

Method Info : Mobile Phase: A: water(0..01%TFA) B: ACN(0..01%TFA)

Gradient: 5%-95% B in 1.5min


Flow Rate: 2.0ml/min

Column: Sunfire C18 3.5um,4.6*50mm


Oven Temperature:50 C

Mass Range: 110-800 UV: 214nm;254nm


*ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.71	0.02	2776138.25	2423237.75	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.70	0.03	2560.68	1310.61	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
1.70	0.02	2333.19	1623.73	100.00

Ret. Time: 1.70

<<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.88

