


HT26769-187A4 DMSO 400.13MHz



Parameter	Value
1 Title	HT26769-187A4
2 Comment	HT26769-187A4
3 Origin	Bruker BioSpin GmbH
4 Owner	nmrsu
5 Site	
6 Spectrometer	spect
7 Author	
8 Solvent	DMSO
9 Temperature	298.4
10 Pulse Sequence	zg30
11 Number of Scans	16
12 Receiver Gain	202
13 Relaxation Delay	1.0000
14 Pulse Width	14.5000
15 Acquisition Time	3.9846
16 Acquisition Date	2012-07-26T14:39:36
17 Modification Date	2012-07-26T14:39:36
18 Spectrometer Frequency	400.13
19 Spectral Width	8223.7
20 Lowest Frequency	-1640.9
21 Nucleus	1H
22 Acquired Size	32768
23 Spectral Size	65536