

File D:\DATA

Acq. Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Analysis Method : C:\Chem32\1\METHODS\1-POS-3MIN.M

Sample Info : *** No target masses specified ***

Method Info : Mobile Phase: A: water(10mM Ammonium hydrogen carbonate) B: ACN


Gradient: 5%-95% B in 1.6min

Flow Rate: 2.0ml/min

Column: XBridge C18,4.6*50mm,3.5um

Oven Temperature: 40 C

MS Range: 110-800



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.44	0.02	290517.00	204199.34	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.45	0.01	1221.95	1597.42	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.45	0.01	1055.19	1396.03	100.00

Ret. Time: 1.44


<<<< POSITIVE SPECTRA >>>>




Ret. Time: 1.98




Ret. Time: 2.00




Ret. Time: 2.04




Ret. Time: 2.06




Ret. Time: 2.10




Ret. Time: 2.20




Ret. Time: 2.21




Ret. Time: 2.26




Ret. Time: 2.28




Ret. Time: 2.32




Ret. Time: 2.33




Ret. Time: 2.36




Ret. Time: 2.38




Ret. Time: 2.40




Ret. Time: 2.41




Ret. Time: 2.43




Ret. Time: 2.45




Ret. Time: 2.48




Ret. Time: 2.50




Ret. Time: 2.51




Ret. Time: 2.53




Ret. Time: 2.57




Ret. Time: 2.60




Ret. Time: 2.62




Ret. Time: 2.64




Ret. Time: 2.71




Ret. Time: 2.73




Ret. Time: 2.75




Ret. Time: 2.82




Ret. Time: 2.84




Ret. Time: 2.92




Ret. Time: 2.97




Ret. Time: 2.99




Ret. Time: 3.07




Ret. Time: 3.15




Ret. Time: 3.19




Ret. Time: 3.22




Ret. Time: 3.29



Ret. Time: 3.32



Ret. Time: 3.39



Ret. Time: 3.42

