

File D:\DATA\T91\MONITOR\NF-080-4-00541-LCMSA033.D

Injection Date : 12 Jun 12 5:17 pm +0800

Tgt Mass (EZX) :

Sample Name : NF-080-4

Location : P1-D-03

Acq. Operator : T91-Monitor

Inj : 1

Spec. Reported : MS Integration

Inj Volume : 2 ul

Acq. Method : D:\METHODS\1-POS-MON.M

Analysis Method : D:\METHODS\1-POS-MON.M

Sample Info : *** No target masses specified ***

Method Info : Column: SunFire C18, 4.6x50mm, 3.5um


Mobile phase: H2O(0.05%TFA) (A) / ACN(0.05%TFA) (B)

Elution program: Gradient from 5 to 95% of B in 1.2min at 2ml/min


Temperature: 50 °C

Detection: UV (214 , 4 nm) and MS (ESI, Pos mode , 70 to 1000 amu)


*ADC1 A, ELSD




*DAD1 A, Sig=214,4 Ref=off



*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Integration Results for ADC1 A, ELSD

RetTim	Width	Area	Height	Area%
1.30	0.03	17041.67	8627.04	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.29	0.03	4104.44	2068.67	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.29	0.03	1399.93	768.61	100.00

Ret. Time: 1.24

<<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.29



Ret. Time: 1.93

