





LC/MS Report

File D:\DATA\MONITOR\211915-06-9-LCMS.D
Injection Date : 22 Oct 12 5:17 pm +0800 Tgt Mass(EZX) :
Sample Name : 211915-06-9 Location : P2-F-05
Acq. Operator : LCMS Inj : 1
Spec. Reported : MS Integration Inj Volume : 2 ul
Acq. Method : C:\Chem32\1\METHODS\1-POS-3MIN.M
Analysis Method : C:\Chem32\1\METHODS\1-POS-3MIN.M
Sample Info : *** No target masses specified ***
Method Info : Mobile Phase: A: water(10mM Ammonium hydrogen carbonate) B: ACN
Gradient: 5%-95% B in 1.5min
Flow Rate: 2.0ml/min
Column: XBridge C18,4.6*50mm,3.5um
Oven Temperature: 40 C
MASS Range:110-800



Integration Results for ELS1 A, Voltage

RetTim	Width	Area	Height	Area%
2.08	0.02	371.76	251.46	100.00

Integration Results for DAD1 A, Sig=214,4 Ref=off


RetTim	Width	Area	Height	Area%
1.61	0.02	16.26	14.21	0.89
1.75	0.01	8.53	10.33	0.46
2.03	0.02	12.28	12.22	0.67
2.07	0.02	1797.11	1761.50	97.98

Integration Results for DAD1 B, Sig=254,4 Ref=off


RetTim	Width	Area	Height	Area%
1.61	0.02	4.43	4.25	0.56
1.75	0.01	4.11	4.89	0.52
2.02	0.02	3.65	3.46	0.47
2.07	0.02	773.18	752.48	98.45

Ret. Time: 1.62


<<<< POSITIVE SPECTRA >>>>



Ret. Time: 1.75



Ret. Time: 1.93



Ret. Time: 2.08

