




*DAD1 A, Sig=220,4 Ref=off




*DAD1 B, Sig=254,4 Ref=off



MSD1 TIC, MS File ES-API, Pos, Scan, Frag: 70



Ion 390, MSD1 390, Target Mass 389 +H Positive, EIC=389.7:390.7



Integration Results for DAD1 A, Sig=220,4 Ref=off


RetTim	Width	Area	Height	Area%
11.31	0.03	2681.45	1417.41	100.00

Integration Results for DAD1 B, Sig=254,4 Ref=off

RetTim	Width	Area	Height	Area%
11.31	0.03	1375.66	735.33	100.00

Ret. Time: 11.31 <<<< POSITIVE SPECTRA >>>>

ES-API Positive





8.405
8.402
8.384
8.382

7.451
7.392
7.388
7.372
7.369
7.261
7.107
7.104

2.938

1.595
1.564
1.549
1.463
1.451
1.379
1.349
1.304
1.287
1.274
1.255
1.237
1.221
0.800
0.781
0.763
-0.000



10 02:15:19 2012